




1 Let's start! Apply conductive gel on both electrodes

C'est parti! Appliquez du gel sur les deux électrodes / Los geht's! Verteile das Gel auf beiden Elektroden / Empecemos! Aplicar gel en los dos electrodos / Vamos começar! Aplique o gel nos dois eletrodos / Cominciamo! Applicare del gel direttamente sugli elettrodi




2 Place leather attachment

Placez l'attache en cuir sur la sangle / Fixiere die Lederhalterung / Coloca el enganche de cuero / Coloque a fixação de couro / Posizionare l'attacco in pelle




3 Place top electrode support under the saddle pad

Placez l'électrode du haut sous le tapis / Platziere die obere Elektrode unter der Schabracke / Coloca el electrodo superior sobre la almohadilla / Coloque o eletrodo superior sob a sela / Posizionare l'elettrodo superiore sotto il sottosella



4 Adjust elastic band & clip-on the snap button

Ajustez l'élastique et clippez le bouton pression / Passe das Elektrodenband an & drücke den Druckknopf zu / Ajusta el elástico y haz clic en el botón de presión / Ajuste a faixa de elástico e pressione o botão de pressão / Aggiustare l'elastico e premere sul bottone a pressione



5 Insert the tracker & you're ready to start!

Insérez le capteur et vous êtes prêt ! / Führe den Sensor ein. Fertig! / Inserta el sensor y ya estas listo para empezar! / Insira o rastreador e você está pronto para começar! / Inserire il sensore per iniziare!

